



J.R. INSTRUMENTS

OLYMPUS

Nikon

Authorized Microscope Service Dealer

PREVENTATIVE SERVICE:

CLEANING: All optical elements / mechanical components will be thoroughly cleaned using manufacturer's recommended procedures.

LUBE: Moving parts requiring lubrication will be degreased with proper solvents and re-lubricated with non-drying lubricants.

OPTICAL SYSTEM: Optical elements will be inspected to insure that all optical alignments (focus adjustment, collimation, resolution & magnification) meet manufacturer's specifications & tolerances.

SYSTEM: An image evaluation will be conducted to insure optimum performance of the microscope and to insure the instrument is optically balanced to full specifications.

LIGHTING: We will align lighting paths, including condensers, Epi-Fluorescence, Bright field Dark field, Phase Contrast, & Polarized illumination. We will also inspect for bad Power cords, lamp housings assemblies & lamps sockets.

WARRANTY: 6 Months labor and workmanship.

J.R. INSTRUMENTS

11838 Jefferson Blvd Culver City, CA. 90230 Tel 310-398-7951
Web site: www.jrinstruments.com service@jrinstruments.com

Fax 310-398-7971
Email: info@jrinstruments.com